SUPAJET

JOINT ASSESSMENT EXAMINATION July/August 2013 **CHEMISTRY DEPARTMENT**

Name	<u> </u>	Class	Class No

233/3 CHEMISTRY PAPER 3 PRACTICAL TIME $-2^{1}/_{4}HRS$.

Index Number:

INSTRUCTIONS.

Answer ALL the questions in the spaces provided in the question paper.

You are NOT allowed to start working with the apparatus for the first 15 minutes of the 21/4 hours allowed for this paper.

This time is to enable you to read the question paper and make sure you have all the chemicals and apparatus that you may need.

Electronic calculators may be used

All working must be clearly shown shown where necessary.

FOR EXAMINER'S USE ONLY.

QUESTION	MAXIMUM SCORE	CANDIDATE'S SCORE
	SCORE	SCORE
1	21	
_		
	11	
	11	
2		
2	0	
3	8	
Total Cooms	40	
Total Score	40	

astRagers.com

- **Q1**. You are provided with:
 - Aqueous sulphuric acid labelled solution A
 - Solution B containing 8.0g per litre of sodium carbonate.
 - Solution C containing sodium hydroxide.

You are required to determine the concentration of solution A, and use the diluted solution to find the concentration of solution C.

Procedure1

Using a picette, place 25.0cm³ of solution A into a 250ml volumetric flask. Add distilled water to make 250cm³ of solution – Label this solution D. Place solution D in a burette. Clean the pipette and use it to place 25.0cm³ of solution B into

conical flask. Add 2 drops of **methyl orange** indicator provided and titrate with solution D. Record your results in the table below. Repeat the titration two more times

and complete the table.

	I	II	III
Final burette reading			
Initial burette reading			
Volume of solution D used (cm ³)			

(4marks)

- a) Evaluate the:
- i) Average volume of solution D used.

(1mark)

ii) Concentration of sodium carbonate in solution B. (R.A.M. Na=23, C=12, O=16)

(1mark)

b) Write the equation of the reaction that occurs (1mark)

.....

- c) Calculate the:
- i) Moles of sodium carbonate reacting with sulphuric acid. (1mark)

	ii) Concentration of sulphuric acid in	n solution D.	(2 marks)	
\$ ⁴ ee	iii) Concentration of sulphuric acid in the	n solution A.		(2 marks)
Cot Wote	Procedure 2.			
	Procedure 2. Place solution D in a burette. Clean solution C into a conical flask. Add and titrate with solution D. Record titration two more times and complete.	2 drops of phen your results in the	olphthalein indi	cator provided
		I	II	III
	Final burette reading			
	Initial burette reading			
	Volume of solution D used (cm ³)			(4marks
	d) Evaluate theAverage volume of s	solution D used.	(1mark)	
	e) Write the equation of the reaction	າ that occurs (1m	ark)	

f) Calculate the :

i) Moles of sodium hydroxide reacting with sulphuric acid. (1mark)

iv) Concentration of sodium hydroxide in solution C in moles per litre.

(2 marks)

Q2. You are provided with:

Solution **F**, which is 1M sodium Hydroxide solution,

Solution **G**, which is 1M hydrochloric acid solution.

You are required to find the Heat of neutralisation of sodium hydroxide.

Procedure:

Measure 50cm³ of solution F into a plastic beaker and record its temperature. Measure 50cm³ of solution G and record its temperature. Pour solution F into the beaker containing solution G. Stir thoroughly and record the highest temperature attained by the mixture. Record your results in the table below.

Temperature of A (°C)	
Temperature of B (°C)	
Highest temperature (°C)	
	(3 marks)

Temperature change ($\Delta\Theta$)	(2 marks)
(Show clearly how you got your answer)	

with liquid Droperty out the tests descri

Q3. You are provided with liquid P; carry out the tests described below.

a) Place a little amount of liquid P in a metallic spatula and ignite it in a Bunsen burner flame.

Observation with	Inferences	
, b		
A'		
age ^{te}		
₹.		
Qasx *		
Ç [©]	(4 ma a rls)	(4 ma a rls)
9×	(1 mark)	(1 mark)

b)Divide the remaining liquid into 3 portions.

i) To the first portion add three drops of acidified potassium dichromate

Observation	Inferences
(1 mark)	(1 mark)

ii) Add an equal volume of distilled water to the the 2nd portion followed by a spatula-end full of sodium carbonate

Observation	·	Inferences	
	(1 mark)		(1 mark)

iii) Add about 3cm³ of ethanoic acid to the 3rd portion followed by 3 drops of sulphuric acid and warm.

Observation		Inferences	
	(1 mark)		(1 mark)