

BIOLOGY PRACTICAL 231/3 2014 MOCKS

CONFIDENTIAL

1. All photographs must be coloured
2. Each candidate to be provided with ;

-3 test tubes

-test tube rack and holder

-a filter paper

-a funnel

-50 ml beaker

Access to

- source of heat

-Benedicts solution

-iodine solution

- Stirring rod

-2 droppers

3. Solution Z prepared by dissolving 50grams of starch and 50grams of glucose in 100mls of warm

distilled water.

Each candidate to be provided with 20 mls