

29.6.3 Chemistry Paper 3 (233/3)

1 You are provided with:

- acid A labelled solution A;
- 2.0 M sodium hydroxide solution labelled solution B;
- Solution C containing 25.0 g per litre of an alkanolic acid.

You are required to:

- (a) prepare a dilute solution of sodium hydroxide, solution B.
- (b) determine the:
 - (i) molar mass of the alkanolic acid
 - (ii) reaction ratio between sodium hydroxide and acid A.

Procedure 1

Using a pipette and a **pipette filler**, place 25.0 cm³ of solution B into a 250.0 ml volumetric flask. Add about 200 cm³ of distilled water. Shake well. Add more distilled water to make up to the mark. Label this solution D. **Retain the remaining solution B for use in procedure II.**

Fill a burette with solution C. Using a clean pipette and a **pipette filler**, place 25.0 cm³ of solution D into a 250 ml conical flask. Add two drops of phenolphthalein indicator and titrate with solution C. Record your results in **table 1**. Repeat the titration two more times and complete the table.

	I	II	III
Final burette reading			
Initial burette reading			
Volume of solution C (cm ³) added			

(4 marks)

Determine the:

- (i) average volume of solution C used; (1 mark)
- (ii) concentration of solution D in moles per litre; (1 mark)
- (iii) concentration of the alkanolic acid in solution C in moles per litre (1 mole of the acid reacts with 3 moles of the base); (1 mark)
- (iv) molar mass of the alkanolic acid. (1 mark)

Procedure II


Fill a **clean** burette with solution A. Place 5 cm³ of solution A into a 100 ml beaker. Measure the initial temperature of solution A in the beaker and record it in **table II**. Using a 10 ml or a 100 ml measuring cylinder, measure 25 cm³ of solution B. Add it to solution A in the beaker and immediately stir the mixture with the thermometer. Record the maximum temperature reached in **table II**. Repeat the experiment with other sets of volumes of solutions A and B and complete the table.

Table II

Volume of solution A (cm ³)	5	9	13	17	21	25
Volume of solution B (cm ³)	25	21	17	13	9	5
Maximum temperature (°C)						
Initial temperature (°C)						
Change in temperature, ΔT						

(6 marks)

- (a) On the grid provided; plot a graph of ΔT (Vertical axis) against the volume of solution A. (3 marks)



- (b) From the graph, determine the volume of solution A which gave the maximum change in temperature. (1 mark)
- (c) Determine the volume of solution B that reacted with the volume of solution A in (b) above. (1 mark)
- (d) Calculate the:
- ratio between the volumes of solutions A and B that neutralised one another; (1 mark)
 - concentration in moles per litre of the acid in solution A. (Assume that the volume ratio is the same as the mole ratio). (1 mark)

2 You are provided with solids E, F and G. Carry out the tests below and write your observations and inferences in the spaces provided.

(a) Place all of solid E in a boiling tube. Add 20 cm³ of distilled water and shake until all the solid dissolves. Label this as solution E.

(i) To about 2 cm³ of solution E in a test-tube, add 4 drops of 2 M sulphuric (VI) acid.

OBSERVATIONS	INFERENCES
(1 mark)	(2 marks)

(ii) To about 2 cm³ of solution E in a test-tube, add 2 M sodium hydroxide dropwise until in excess.

OBSERVATIONS	INFERENCES
(1 mark)	(1 mark)

(iii) Place one half of solid F in a test-tube. Add 2 cm³ of distilled water and shake well. Add 4 drops of this solution to about 2 cm³ of solution E in a test-tube.

OBSERVATIONS	INFERENCES
(1 mark)	(1 mark)

(iv) To about 2 cm³ of solution E in a test tube, add 2 drops of aqueous potassium iodide.

OBSERVATIONS	INFERENCES
(1 mark)	(1 mark)

(b) (i) Using a metallic spatula, ignite about one half of solid G in a Bunsen burner flame.

OBSERVATIONS	INFERENCES
(1 mark)	(1 mark)

(ii) Place the other half of solid G into a boiling tube. Add 15 cm³ of distilled water and shake well. Label this solution G. Use this solution for the following tests.

I Place 2 cm³ of solution G in a test-tube and determine its pH.

OBSERVATIONS**INFERENCE**

(1 mark)

(1 mark)

II To about 2 cm³ of the solution obtained in (ii) above, add 3 drops of acidified potassium manganate (VII).

OBSERVATIONS**INFERENCE**

(1 mark)

(1 mark)

III To about 2 cm³ of the solution obtained in (ii) above, add 2 drops of bromine water.

OBSERVATIONS**INFERENCE**

(1 mark)

(1 mark)

(iii) To the remaining solution G in the boiling tube, add the other half of solid F.

OBSERVATIONS**INFERENCE**

(1 mark)

(1 mark)